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Key  
O Number of ring marks on the attenuation element  
Mf: Mean Attenuation in dB  
Sf: Standard Deviation in dB  
APV: Assumed Protection in dB  
SNR: Single Number Rating in dB  
L: Attenuation values for low frequency noises in dB  
M: Attenuation values for medium frequency noises in dB  
H: Attenuation values for high frequency noises in dB  
mean: SNRm, Hm, Mm, Lm  
sd: SNRs, Hs, Ms, Ls  
Result: Result

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GB

#### User instructions for hearing protection SILOC

With your new, custom-made hearing protection SILOC you have decided to protect your hearing against damaging or disturbing noise by using custom-made protectors. Your hearing protection device has been manufactured to fit your ear perfectly and is guaranteed to offer you the highest possible wearing comfort. The hearing protector is worn in the ear. It can be inserted with minimal effort, is highly compatible with other protective equipment (e.g. protective glasses) and above all, is appropriate for long stays in noise zones. Properly used, it protects you from hazardous noise. No other use is approved and is therefore not allowed. We hope you enjoy it thoroughly!

#### Prior to using your hearing protection device

Please read these instructions carefully to familiarize yourself with the proper handling of the hearing protection device. Only if the following suggestions are adhered to can effective protection and the durability of your hearing protectors be guaranteed. Keep these instructions in a safe place.

#### Standard features of your hearing protection device

- Hearing protection earmold (Note color marks: right ear/ red, left ear /blue)
- Attenuation element for the hearing protectors (1 pair)
- Wax pin (use the metal loop to remove cerumen deposits from the hearing protectors)
- Case
- User instructions and warranty card

#### Additional features of your hearing protection device

- Connecting cord (detachable connection between the right and left hearing protector)
- Safety clip to attach connecting cord to clothes
- otosense care gel to ease inserting the protectors
- Antimicrobial coating

#### Optional additional features and products

- Attenuation elements with different attenuation effect – “Technical data” (see table)
- care products for cleaning and disinfecting (e.g. otosense disinfectant spray)

Your hearing protection device is equipped with a sound canal, which contains an attenuation element. Depending on your requirements, you may use other attenuation elements, which can be purchased from your specialist.

#### Inserting and removing your hearing protection device (see figure 2 to 4)

- With clean hands, take the hearing protection SILOC out of the packaging.
- Make sure that your entire hearing protectors have not been damaged and are clean. Do not insert a damaged hearing protector in your ear; in case of damage, please contact your specialist.
- Note the red and blue color marks on your hearing protectors: red = protector for the right ear; blue = protector for the left ear.
- Take the right protector into your right hand, the left protector into your left hand. Be sure to grasp the protector from the bottom with your thumb and from the top with your index finger (figure 2).
- Take the hearing protector to your ear without twisting it. Angling it from behind, insert it into the auditory canal (figure 3). Be sure to insert the hearing protection device deeply enough and position it comfortably (figure 4).
- To facilitate the insertion of the protectors, we recommend using otosense care gel (see additional products).
- To remove the hearing protector grasp the handle with your thumb and index finger

and pull the protector out of your ear. Alternatively you can free the hearing protector by lifting it slightly with your thumb and then removing it from your ear with your thumb and index finger (see figure 2).

#### WARNING

##### A) Noise impact on hearing without hearing protection

Even a brief exposure to a high noise level while wearing improperly inserted hearing protectors or no hearing protection devices at all may damage your hearing.

##### B) Protective effect that is too low when the hearing protector is incorrectly inserted in the ear

Only expose yourself to high noise levels while wearing your custom-made hearing protectors after you are familiar with handling the protectors and inserting them properly into the auditory canal. Have the correct use shown to you e.g. by the manufacturer or the supplier of the hearing protection earmold or your physician. It may be helpful for beginners to practise the correct technique in front of a mirror.

##### C) Handling and daily checks

The sound canal must be sealed using the original stainless steel attenuation element, and may only be removed by a specialist for the purposes of cleaning or testing. If the stainless steel attenuation element is missing, do not wear the hearing protection devices in extreme noise levels! Please check your hearing protector for this on a daily basis.

##### D) Leaking /reduced protective effect

Anatomical changes in the shape of the ear may lead to the deterioration in secure earmold fit, e.g. when there is a change in weight or there have been ear operations. Safe protection is only achieved when a specialist has tested the hearing protection earmold after it is delivered and then tests it again at regular intervals of at most three years.

##### E) Reduced perception

Please consider: While wearing the hearing protection devices, the signals in the

surrounding area will also be muffled (e.g. warning signals). In addition, be sure to adhere to the regulations regarding the use of hearing protection devices in road traffic/track superstructure/railway operation.

##### F) Risk of injury with improper handling

Do not pull the hearing protectors out of the ear by the connecting cord. Sudden or fast removal of the ear plugs out of the auditory canal may damage the ear drum.

##### G) Risk of the hearing protector being pulled out.

These ear plugs should not be used where there is a risk that the connecting cord could be caught up during use. We recommend detaching the connecting cord with the connecting fasteners from the hearing protection device. After the connecting cord has been detached, it can only be reattached by inserting it again (note the connecting fasteners in figure 5).

##### H) Connecting additional equipment

As a hearing protection device the SILOC may only be used in connection with additional products (e.g. attenuation element) that have been type-tested with SILOC. More information is available from your specialist or the manufacturer.

##### I) Skin reaction

The hearing protector is exclusively manufactured from materials that have been proven over the years, tested and that are bio-compatible. If you were nonetheless to notice symptoms, remove the hearing protector immediately and inform your doctor, specialist consultant or supervisors.

##### J) For small children there is a swallowing hazard!

The product must be kept out of reach of small children.

##### K) Ageing and wear

Custom made earmolds and attenuation elements are designed for daily use. However, we recommend having a new spare earmold made after at most 6 years. If necessary, the attenuation element may continue to be used after successful funtion testing.

#### Storage and care

This product may suffer if certain chemical substances are used. Your hearing protectors must be cleaned regularly to avoid a decrease in its protective effectiveness, or the development of skin irritations or other ear problems. Use either special cleaning products offered by the manufacturer or at least water and soap. Dry the hearing protectors thoroughly before inserting them again. Do not use chemical solvents, abrasive cleansers or brushes. The attenuation elements must be cleaned professionally. More information is available from your hearing protection professional or the manufacturer. Store the SILOC in the case which is included in delivery.

#### ATTENTION

These products may not be altered or modified. If the suggestions in these instructions are not adhered to, the effectiveness of the hearing protection device may be affected considerably. Only by handling the hearing protectors carefully and wearing the hearing protection SILOC continuously in extreme noise levels, are you protecting your hearing as best as possible. Therefore, we urgently recommend using and maintaining the hearing protection SILOC according to these suggestions.

#### CE Declaration of Conformity

SILOC hearing protection devices are manufactured individually after ear impression taking. Certified and monitored (Module D) by PZT GmbH, Bismarckstrasse 264 B, D-26389 Wilhelmshaven. European notified body with the identification number 1974. otofaktor GmbH, declares herewith that the Personal Protective Equipment SILOC corresponds to the PPE Regulation (EU) 2016/425. The complete text of the EC Declaration of Conformity is available at the following Internet address: <https://www.otofaktur.de/files/otf/documents/EU-Konformitaetserklaerung.pdf>

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